

XQ-320

Excellence in quality and usability





Lighting the way
with diagnostics

SWCA will be the reliable diagnostics company in West and Central Africa providing medical laboratories with quality testing, operational efficiency, and access to innovative diagnostics that ensure public safety and optimal output from hospital systems.

www.sysmex-wca.com

01

XQ-320 overview and specifications

Concept

Convenient



Faster



Slimmer



Lighter



XQ-320: Fusion of requirements and inheritance

inherit

- XN-Series: Data correlation
- XN-Series: High-end quality
- XN-L Series: Software UI
- XN-L Series: Hardware
- XP-300: Technology, parameters
- XP-300: Reagents



evolve

Market requirements

- QC lockout
- Flexible user management
- Barcode login
- Assay sheet 2D barcode
- Understandable flagging names
- Higher throughput
- Lower consumption

- XQ-320 represents the 'best' of Sysmex analysers to answer even better to your needs

Laboratory management

QC lockout

Host format XQ/HL7

IPU screen lock (≥ 1 min)

Flexible user management

Up to 500 users

Barcode logon

QC measurement permission

Ease of use

5 USB ports (1 on front)

Wi-Fi/Bluetooth connectivity

Assay sheet 2D barcode

Faster wakeup

Additional languages

Histograms next to numerical results

Understandable flagging names

Less user maintenance

Specifications

Higher durability

Higher throughput

Lower aspiration volume

Smaller bench size

Larger temperature range


Larger humidity range

Larger memory

Larger touch screen

Lower reagent consumption

Lower electricity consumption



XQ-320 at a glance

Integrated IPU

Colour LCD with touch screen

Reagent holder

Power switch



Analysis status indicator

Sample start switch

XQ-320: Built for convenience



Convenient

- 10.4" display
- Numerous new software features
- Integrated IFU



XQ-320: Built for fast results




result in
<60
seconds



Faster

- Short time to result
- Decreased TAT



XQ-320: Compact and robust hardware



XQ-320: Compact and robust hardware



Lighter

↓
22 kg



↓
30 kg



Specifications in comparison to XP-300: Overview



Improved usability

	XP-300	XQ-320
Aspiration volume WB mode	50 μ L	16 μL
Aspiration volume PD mode	200 μ L	65 μL (diluted blood)
Ambient temperature	15–30°C	15–30°C
Humidity	30–85%	20–85%
Interface [USB LAN Serial Parallel]	- 1 1 1	5 2 1 -
Electric power	< 200 VA	< 160 VA
Network	Wired	Wired Wireless

Specifications in comparison to XP-300: Memory



Enhanced memory function

	XP-300	XQ-320
Sample storage	40,000	500,000
Patient information	Sample ID only	10,000
Registered wards	—	200
Registered doctors	—	200
QC files	6 files 60 plots/file	99 files 300 plots/file
Calibration log	—	20
User management	Operator ID only	500 users*

* Including customisable settings and barcode identification

Specifications in comparison to XP-300: Maintenance




Less user maintenance, increased uptime


	XP-300	XQ-320
Startup	Startup	Startup, auto wakeup
Shutdown	Shutdown + Cellclean	Power off, sleep mode
Cleaning of SRV tray	Weekly	—
Cleaning of SRV	Every 3 months or 4,500 tests	—
TD cleaning	Monthly or 1,500 tests	—
Waste chamber cleaning	Monthly or 1,500 tests	When needed
Routine cleaning	—	Weekly or 400 tests

Power on/Wakeup

- Sequences at startup or wakeup are much faster
- XQ-320 is quickly ready for measurement
- First Sysmex haematology analyser with automatic wakeup schedule function



New features: Connectivity




- Connectivity via ethernet, Bluetooth or Wi-Fi
- HL7 protocol ensures seamless integration into hospital networks



Result screen


Everything at a glance

- Parameters, histograms and messages on one screen
- Easy to read flags
- 20 diagnostic parameters



Diagnostic parameters

Channel	Parameters
WBC*	WBC, NEUT#/%, LYMPH#/%, MXD#/%
RBC/PLT	RBC, HCT, MCV, MCH, MCHC, PLT, RDW-SD, RDW-CV, PDW, MPV, P-LCR, PCT
HGB	HGB



* WBC parameter names can be adapted: LYM% (W-SCR), MXD% (W-MCR), NEUT% (W-LCR), LYM# (W-SCC), MXD# (W-MCC), NEUT# (W-LCC)

Reagents




Function	Reagent name	Standard	Option
Diluent	Cellpack	1 x 20 L (politainer)	1 x 10 L (politainer)
Lysing reagent	Stromatolyser-WH	3 x 500 mL (bottle)	
Detergent	Cellclean	Cellclean (1 x 50 mL)	

Caresphere™ XQC - interlaboratory comparison

Caresphere™ XQC

- Caresphere™ XQC judges the performance of an individual analyser by comparing its QC results to peer group data.
- Outliers are quickly identified.
- No additional measurements or materials are needed.



Summary of XQ-320 benefits

Some Benefits of the XQ-320 over XP-300 at a Quick Glance

	XQ-320	XP-300
Reagent consumption	Consumes relatively less Cellpack=30ml Stromatolyser=0.8ml Cell clean= once/week	Consumes relatively more Cellpack=34ml Stromatolyser=1ml Cell clean= every shutdown
Throughput	Higher through put -Relatively Faster Whole blood mode=70 tests/hour	Lower through put -Relatively Slower Whole blood mode=60 tests/hour
Maintenance	Less maintenance activity (apart from start-up and shutdown) – routine cleaning	More maintenance activity (apart from start-up and shutdown) - Cleaning of SRV tray; Cleaning of SRV; TD cleaning
Sample Aspiration volumes	Requires relatively less sample volume: Whole blood - 15ul Predilution – 65ul	Requires relatively more sample volume Whole blood - 50ul Predilution – 200ul
Weight	Relatively lighter 22kg	Relatively heavier 30kg
Dimensions	Relatively slimmer 365(W) x 440(D) x 450(H) mm	Relatively larger 420(W) x 355(D) x 480(H) mm
Data Storage	Stores more data: Sample data - 500, 000 samples including histograms QC files - 99 files 300 plots/file Patient info - 10,000	Stores relatively less data: Samples data-40,000 samples including histograms QC files - 6 files 60 plots/file Patient info - Sample ID only
Network connection	Wired and Wireless	Wired


02

Market Analysis and product strategy

3 PD market in West Africa

Market/ Environment Analysis	<ul style="list-style-type: none">• Biggest than 5 PD market in volume• Primary healthcare - Medium and low size hospitals and private clinics• Remote country• General medicine (+++), less specialities• 15- 30 samples /day on average• Price sensitive market (+++)
Competitor Analysis	<ul style="list-style-type: none">• Lot of Chinese brands: Mindray 3PD, Rayto 3PD, URIT 3PD• Manual count
Stakeholder Analysis	<ul style="list-style-type: none">• Hospital management: Medical doctors or pharmacists are decision makers (+++)• Lab manager: generally involved in the purchase process as advisor or decision makers• Others: procurement managers, director...• MOH in case of tenders

Behaviour change




Product strategy

Vision	Best in class in the 3 PD market in West Africa to answer more than customers needs			
Positioning	3PD – part differential with a dedicated neutrophil count for primary healthcare in Africa			
Strategic Imperatives	Create awareness about burden of infectious disease	Distributors sales and app teams trained	Replacement strategy / tenders	Product accessibility
Required stakeholder behavior change	<p>Burden and management of infectious disease (LRTIs in children)</p> <p>By GPs and clinicians CMEs in primary care settings(training).</p>	<p>Distributor teams trained about our key differences vs competitions</p> <p>To replace competition instruments in the lab</p> <p>By .showing the added value of our 3 PD analysers</p>	<p>Product features: low maintenance (+++), low aspiration volume, caresphere</p> <p>By selling our service support (technical, application and scientific)</p>	<ul style="list-style-type: none"> • Same price than XP-300
SUCCESS METRIC	<ul style="list-style-type: none"> • GPs influence lab management to purchase . • Reagent consumption by test prescription 	<ul style="list-style-type: none"> • Market share gained from competition > 10 % first year. 	Perception of Sysmex as a MOH partner of choice	Number of units sold / year

Commercialization strategy

key message



Application and scientific Affairs development strategy

key message

Key challenges

GPs and clinicians to interpret FBC parameters

Unmet needs / Gaps

Lack of knowledge about diagnostic parameters that can help them better treat patient with infections

Lack of application and clinical training .

Strategy

- .Communication and CMEs campaign on dedicated neutrophils count.
- Distributors to hire dedicated application staff to cover customers and clinicians training in-country

Summary XQ-320 WCA Strategy


- **WHERE-** Primary Healthcare Centres (MOH & Privates) / Satellite labs
- **HOW;**
 - Local strategies to gain market share in the 3 PD market
 - » Clinicians campaign focused on the importance of dedicated neutrophils count in the management of infection in primary care
 - » End user campaign: XQ features, ease of Use, low sample volume and maintenance, Large screen
 - Marketing launch activities in Nigeria, IC, Senegal, Benin DRC, Congo , Gabon , & Cameroun

03

SPMs

GPs and clinicians campaign: How to differentiate infection with a dedicated neutrophils count

XQ-Series



Primary healthcare diseases in Africa

In sub-Saharan Africa, patients with infectious diseases are frequently encountered in general medicine and paediatrics departments. Lower respiratory tract infections (LRTIs) represent a significant cause of mortality, ranked third after HIV/AIDS and malaria in sub-Saharan Africa overall, and first in nine African countries [1].

LRTIs are commonly caused by bacteria, including *Streptococcus pneumoniae* (*pneumococcus*), *Haemophilus influenzae*, *Klebsiella pneumoniae*, *Staphylococcus aureus*, *Aerobacter* species, *Streptococcus mitis*, *Parabacterium aeruginosus*, *Escherichia coli* and *Proteus* species. LRTIs include bronchitis, pneumonia, and tracheolitis.

The threat posed by pneumonia

Pneumonia is a form of acute respiratory tract infection that affects the lungs. The lungs are made up of small sacs called 'alveoli', which fill with air when a healthy person breathes. When an individual has pneumonia, the alveoli are filled with fluid, which makes breathing painful and limits oxygen intake.

Pneumonia is the single largest infectious cause of death in children worldwide. Pneumonia killed 740 180 children under the age of five in 2019, accounting for 22% of all deaths in children aged one to five. Pneumonia affects children and families everywhere, but deaths rates are highest in South Asia and sub-Saharan Africa. It can be caused by viruses, bacteria, or fungi. Children can be protected from pneumonia, as it can be prevented with simple interventions, and treated at low cost, if it is properly diagnosed [2].

Pneumonia caused by bacteria can be treated with antibiotics, but only one third of children with pneumonia receive the antibiotics they need.

References

- [1] A large-scale seroprevalence study of respiratory tract infections in people with low average literacy rates in Africa: the African seroprevalence of disease and health (ASD) study. *PLoS One*. 2019; 14(12): e0224000.
- [2] Pneumonia and child health: <https://www.who.int/news-room/fact-sheets/detail/pneumonia>

1 | XQ-Series

How do we differentiate infections?

Clinical symptoms



Laboratory investigations



In a patient with an acute febrile illness, infection must always be considered. The CBC including a dedicated neutrophil count is helpful in guiding the diagnosis towards a bacterial or viral infection, providing the result in only one minute. With this information, physicians will be better able to decide whether or not they need to prescribe any antibiotics.

Over low of CBC parameter result ranges characteristic for viral and bacterial infections:

Parameter	Reference range	Viral Infections	Bacterial Infections
White blood cells	4-10 × 10 ⁹ /L	>10 × 10 ⁹ /L	>10 × 10 ⁹ /L
Neutrophils	16-7 × 10 ⁹ /L	≤ 7 × 10 ⁹ /L	> 7 × 10 ⁹ /L
Lymphocytes	15-4 × 10 ⁹ /L	> 4 × 10 ⁹ /L	≤ 4 × 10 ⁹ /L

Note: Laboratory diagnostic results must be reviewed by a physician, considering the clinical case and age, and the clinical concerns. As with any case, other laboratory tests and examinations are also necessary to confirm the diagnosis of an infectious disease. Please refer to the relevant literature for detailed information on viral and bacterial infections and related laboratory tests and parameters.



4 | XQ-Series

Clinical case: Meet Kwasi

A mother came to the district hospital with her 3-year-old son Kwasi. Upon presentation to the physician, she complained about his cough. The physical examination showed a raised body temperature of 39.5°C and a weight of 10 kg. When examining the lungs, the physician noted rales during breathing, and tachypnea. A blood sample was drawn and analysed.

Laboratory diagnostics

The CBC investigations revealed the following results.

CBC parameters	Kwasi's results	Reference ranges
Haemoglobin	12 g/dL (normal)	12 - 16 g/dL*
Neutrophils	18 × 10 ⁹ /L (high)	1.6 - 7 × 10 ⁹ /L*
Lymphocytes	1 × 10 ⁹ /L (low)	1.5 - 4 × 10 ⁹ /L

*Reference ranges for children



The white blood cell results with their high neutrophil and low lymphocyte counts pointed to a bacterial infection. After further investigations, such as testing for other biomarkers like CRP and obtaining a chest radiograph, it was confirmed that Kwasi had a community-acquired pneumonia. Treatment could then begin without further delay.



1 | XQ-Series

XQ-Series

We have listened

XQ-320 represents the best of other Sysmex analysers to answer even better to your needs. While keeping the proven technology of the XP-320, we have added new features used in our recent 5-part differential analysers, such as a large touch screen, the user interface and the flagging messages. What's more, results correlate well also with the flagship XN-Series, making the XQ a match in any lab network. From sample aspiration to result presentation, XQ supports labs with its convincing features and benefits.

- ✓ The XQ-320's robustness, combined with its low demand for maintenance, results in high instrument uptime.
- ✓ The low sample volume of just 30 µL makes challenging samples, e.g. from elderly patients or children, much easier to process.
- ✓ Advanced QC settings support you in fulfilling your regulatory requirements.
- ✓ The various possibilities for connecting XQ-320, e.g. to your network, the Extended IPI, or CareSphere™ – and the close correlation with Sysmex 5-part differential results – make XQ-320 a perfect partner in satellite laboratories within a hospital or lab network.



Worry less, gain more

Technology should support you reliably in what matters most: caring for your patients. It should not require much attention or leave results open for doubt. That's why we took care of the easiest possible operation and result presentation, combined with compact and robust hardware, all backed up by our renowned service and training.

- ✓ Space for medical devices is often limited. We understand that – the space-conscious design of our XQ-Series takes this into account. The slim footprint leaves you more room on your lab bench.
- ✓ Your patients are eagerly waiting for their results – with XQ-Series, they are available in less than an minute, achieving a throughput of up to 70 samples/hour.
- ✓ The clear flagging messages provided by XQ-320 support you in interpreting the results.
- ✓ Several types of printers can be connected via USB, Wi-Fi or Bluetooth, allowing flexibility.



Quick access, safe technology

Time is of the essence and regulations are becoming more and more strict, especially in working areas like ERs, ICUs or other rapid-response scenarios. XQ finds the balance between these requirements that make those areas special. With many different users operating the systems, we've made sure login is simple without compromising on safety. And results are presented at a glance with clear information, so you know quickly where the patient's at.

- ✓ Quick access to results is especially important in an emergency. With easy login via barcode and the quick start-up function, your XQ analyser is ready in an instant.
- ✓ Numerical results, histograms and flag messages are presented together on one screen, providing the complete result at a glance.
- ✓ The new XQ-320 features support you in fulfilling the increasing regulatory requirements, e.g. for accreditation. Automatic logout ensures patient data safety in easily accessible areas, such as ERs. The QC lockout prevents producing results that do not comply with regulatory QC requirements.
- ✓ The CareSphere™ software for the monitoring of analyser and QC status and the optional HLT protocol ensure seamless integration of XQ analysers into your hospital network.



5 USB ports
2 LAN ports

10.4" touchscreen

Scan all assay values in one go

QC lockout function

Connectivity via Wi-Fi or Bluetooth

Easy-to-read stand flags

Results in less than 60 seconds

Complete result on one screen

Screen lockout function

Easy login by scanning barcode/QR code

HL7 protocol for data exchange

04


Sysmex academy training

XQ-320: Training

- **For everyone including Sales colleagues**
 - XQ-320 Basic Online Training on Sysmex Academy platform <https://wca.sysmex-academy.com> (for every one including potential clients)
 - Video Tutorial: XQ-320 basic training to register via LMS.
 - Haematology Sales and Marketing training will be offered to Sales representatives in our academy.
- **For Field Service Engineer**
 - Additional to XQ-320 Online Training and Video Tutorial: XQ-320 basic training Sysmex Academy platform
 - ‘XQ Technical Introduction’ online training on LMS for all Engineer already trained SGST 2
 - New technical colleagues will be trained during SGST 2
- **For Application Colleagues**
 - Additional to XQ-320 Online Training and Video Tutorial: XQ-320 basic training
 - XQ Advanced Application Training. Will be available in LMS

XQ-320: More than just an analyser!

Knowledge
& training



Scientific &
Application
support



Technical
support



Thank you for your attention!

XQ-320



Lighting the way **with diagnostics**

